


Fitting Guide

Nathos Nova M Slimtube

Nathos™ Nova M Slimtube


- 1 Right or left indicator* on back of hearing aid
- 2 Battery compartment with ON/OFF switch
- 3 Multi-function button
- 4 Slim tube
- 5 Retention cord
- 6 Speaker (with Dome earpiece attached)




*Right = red | Left = blue

Turning on the hearing system and inserting into your ear

Switch on/off the hearing system by opening/closing the battery compartment (below).




on




off


Changing the battery



Remove the sticker from the new battery and wait two minutes.




Open the battery door.




Place battery in the battery door with the "+" symbol facing upwards.

Inserting the hearing aid



Place the hearing aid behind your ear.



Insert the earpiece into your ear canal.



If there is a retention cord attached to the earpiece, tuck it into the bowl of your ear to secure your hearing aid.

Changing the volume/program

The Multi-function button functions as a volume control and/or a program change depending on the hearing aid programming.

